

IQ+FLOW

IQPD-600C

Microfluidic Forward Pressure Controller, Downported



Technical specifications

Measurement & control

Type of media	Chip sensor: dry, clean, non-flammable and non-corrosive gases. Media isolated sensor: gases compatible with SS316L and Viton
Accuracy	$< \pm 0.5\% \text{ FS}$
Repeatability	$< \pm 0.1\% \text{ FS}$
Response time (sensor)	$\approx 95\% \text{ 5 msec}$
Operating temperature	$5 \dots 50^\circ\text{C}$
Temperature sensitivity	span: $< 0.1\% \text{ Rd}/^\circ\text{C}$; zero: $< 0.05\% \text{ FS}/^\circ\text{C}$
Leak integrity, outboard	$1 \times 10^{-6} \text{ mbar-l/sec He}$
Pressure range sensor	min. 0.01 - 0.5 bar, max. 0.2 - 10 bar
Max. Kv-value	2.37×10^{-3}

Approvals

Marking	CE
---------	----

Mechanical specs

Pressure rating (PN)	0
Ingress protection	IP40
Material wetted parts	body: aluminium or stainless steel SS316L; chip sensor (flow or pressure): Si, SiOx, epoxy, aluminium; media-isolated sensor (pressure): SS316L
Sealing material	FKM/Viton® seals and plunger (default); FKM/Viton® seals and FFKM/Kalrez® plunger (option)
Process connections	downported construction
Weight	120 g (Aluminum) / 180 g (SS316L)

Electrical properties

Power supply	+15...+24 Vdc, 50 mA for meter; add 50 mA for control valve
Analog output	0...5 (10) Vdc or 0 (4)...20 mA (sourcing)
Analog setpoint	0...5 (10) Vdc or 0 (4)...20 mA (sinking)
Digital communication	RS232/RS485 (Modbus RTU/ASCII or FLOW-BUS)

Electrical interfaces

Power (main connector)	RJ45
Function (main connector)	RS232, RS485, Analog
Modbus RTU/ASCII/FLOW-BUS	RJ45